

# CelCradle-500 Technical Report XVII

## Cultivation of Sf-21 Insect cells in TNM-HF/5%FBS media

### Table of Contents

1 Description .....	2
2 Material .....	2
3 Protocol .....	2
3.1 Inoculum Preparation .....	2
3.2 Preparation before cell seeding.....	2
3.3 Inoculation.....	2
3.4 Culture .....	3
4 Result.....	3
5 Summary .....	5
6 VacciXcell Technical Support .....	6

## 1 Description

CelCradle-500 provides a powerful tool to achieve high cell density and high productivity of target bioproducts in a cell culture because it has a unique feature of offering high oxygen transfer and low shear stress culture environment. Users can easily collect highly concentrated cells, virus or secreted products from one 500 ml CelCradle-500 bottle. In this study, the application of CelCradle-500 culture medium is illustrated.  $2.5 \times 10^8$  for growth of Sf-21 insect cells in TNM-FH/5%FBS SF-21 Insect cells were seeded in one CelCradle-500 unit. A total of 16.5 folds increase and  $4.3 \times 10^9$  cells can be reached at 9 days culture by consuming 1.5 L culture medium. This technical sheet provides a general protocol for users to start up their culture. However, the optimum condition of each cell culture for each case may require the users to determine.

## 2 Material

Device	Cell Line	Medium	Seed
CelCradle-500	Sf-21	TNM-FH/5%FBS	$2.5 \times 10^8$ cells/bottle

## 3 Protocol

### 3.1 Inoculum Preparation

Prepare one 250 ml spinner flask and inoculate  $3.0 \times 10^5$  suspend cells/ml in 120 ml TNM-HF/10%FBS culture media. Culture at 75 rpm, 28°C for 4 days. After cell density reaches above  $1.5 \times 10^6$  cells/ml and viability remain above 95%, it is ready for the preparation of inoculation. Collect  $2.5 \times 10^8$  suspend cells from the spinner flask by centrifugation and collect in one 50 ml centrifuge tube with 30 ml fresh media.

### 3.2 Preparation before cell seeding

Place CelCradle Stage controller in a 28°C incubator. Set up the inoculation parameters (See below). Warm up TNM-HF/5%FBS medium in 28°C water bath. Take out one CelCradle bottle aseptically and place it in a biosafety cabinet. Open the cap and add 470 ml fresh culture medium in the bottle.

### 3.3 Inoculation

Open the cap and distribute 30 ml media containing  $2.5 \times 10^8$  suspended cells that has been prepared previously on top of the matrix of CelCradle-500. Bring the bottle and lock up on the CelCradle Stage controller in incubator at 28°C and start the run immediately. Avoid swirling or shaking the bottle before compression.

### 3.4 Culture

Press “START” button to start the controller. After 3 hours, reset the parameters for culture condition. Usually, above 90% cells will be immobilized in the matrices within 30minutes. The inoculation parameters are set as below:

Rising rate	Top Holding Time	Down Rate	Bottom Holding Time
2.0 mm/s	20 sec	2.0 mm/s	0 sec


The culture parameters are set as below:

Rising rate	Top Holding Time	Down Rate	Bottom Holding Time
2.0 mm/s	0 sec	2.0 mm/s	1 min 30 sec

Due to there are only 0.7 g/L glucose in the culture medium, the pH is very stable during whole culture period. Culture medium were replenished every 3 days and carrier samples were took to count and observe the cell growth.

## 4 Result

### Cell Growth Profile




The cell density started from  $2.5 \times 10^8$  and reached  $4.3 \times 10^9$  with 9 days culture. The growth rate start to slow down after the 6<sup>th</sup> day, which means the medium exchange frequency should be increased from every 3 days to every 2 days.

### Cell Count on carriers


Table 1. The cell number on carriers (CVD method)

Cultivation Time		CVD Nuclear Count		
Days	Hours	2-3 samples raw data	Average	SE
0	0	*	$2.5 \times 10^8$	0
3	76.5	$8.65 \times 10^8$	$1.06 \times 10^9$	$1.93 \times 10^8$
		$1.25 \times 10^9$		
6	148	$2.80 \times 10^9$	$2.80 \times 10^9$	0
		$2.80 \times 10^9$		
8	217	$4.08 \times 10^9$	$3.71 \times 10^9$	$2.03 \times 10^8$
		$3.38 \times 10^9$		
		$3.68 \times 10^9$		


*pH*




pH is stable and even raise during culture. Maintenance of culture is straightforward.



Cells at Day 3



Cells at Day 6



Cells at Day 9

The result indicates that CelCradle cell culture system can be applied in Sf-21 insect cell culture for high cell density culture with very little care on pH control and nutrient supply.

## 5 Summary

Seed	Inoculum Volume	Medium Volume	Medium
2.5 x 10 <sup>8</sup> cells/bottle	30 ml	500 ml	TNF-FH/5%FBS
Total Culture Age	Total Medium Consumed	Total Medium Replenish Frequency	Max. Cell Density
217 hours	1.5L	2 times	4.33 x 10 <sup>9</sup> cells/bottle

## 6 Vaccixcell Technical Support

For queries and comments, please contact the Vaccixcell Technical Support team.

Email: [mail@vaccixcell.com](mailto:mail@vaccixcell.com)

Address: 21 Changi South Street 1, Singapore 486 777

Telephone: +65 6542 0833

Website: <http://vaccixcell.com/>