

**ESCO**  
HEALTHCARE



# MiniTide<sup>®</sup>

The Tiny Titan of  
Adherent Bioprocessing



# MiniTide®

## The Tiny Titan of Adherent Bioprocessing

MiniTide® bioreactor leverages the proprietary Tide Motion® technology in a compact and user-friendly design. This facilitates linear scalability for research-oriented experiments. The streamlined bioprocess design expedites licensing due to the consistent application of the core Tide Motion® technology across all scales. This consistency ensures a seamless transition from initial discovery stages to in-vivo studies, GLP compliance testing, clinical trials, and ultimately, commercialization. This approach optimizes the patent lifespan by employing the same core technology throughout the entire development pipeline.

Esco MiniTide® simplifies the high-throughput screening of exosomes derived from diverse biological sources, covering primary tissues, cancer cells, and immortalized cell lines. This capacity is especially advantageous for the development of diagnostic and screening assays. It acknowledges the substantial variability in therapeutic potential between exosome populations, potentially influenced by culture conditions. The capability of MiniTide® to generate exosomes through cell culture enables researchers to investigate these contributing factors and pinpoint functionally potent exosome subsets.

**Tide Motion Platform**  
(Schematic with Cell Harvesting System)

**CelCradle X® Automated Cell Harvesting System (CCX-ACHS)**  
(CelCradles X Bottles)

**TideXcell® Automated Cell Harvesting System (TXC-ACHS)**

- Closed harvesting system for:
  - ✓ Biomass Expansion
  - ✓ Autologous/Allogeneic Cell Therapy
  - ✓ Intracellular Virus

MiniTide®      CelCradle X®      TideCell / TideXcell® 2 L      TideCell / TideXcell® 20 L      TideCell / TideXcell® 100-300 L      TideCell / TideXcell®

Tide 2 mm/s: 3,000-5,000 L (stainless steel tanks)  
Tide 6 mm/s: 55,000 L

3D seed to 3D bioreactor unlike many other bioreactors which seed in 2D

- GMP (TideCell®) and cGMP units (TideXcell®)
- Surface area:  $3 \times 10^5 \text{ cm}^2$  (110 grams of macrocarriers)
- Larger scale manufacturing


### Features

- Utilizes proprietary Tide Motion® technology in a compact, user-friendly design
- Disposable, single use vessel
- Siemens HMI PLC-based monitoring system
- Enables efficient research with limited samples, costly media, and expensive analyses
- Facilitates rapid exosome generation from primary tissues, cancer cells, and immortalized cell lines

### Applications

- Supports exosome-based research for therapeutic and diagnostic development
- Culture of anchorage-dependent cells with BioNOC II® and BioMESH® macrocarriers
- Proof-of-concept for 3D carrier culture
- MSC/EVs study or process optimization
- Vaccine/virus study or process optimization
- Tissue engineering

## The Motion Principle of MiniTide®



MiniTide® functions based on the Tide Motion principle, where cells adhere to macroporous carriers and are alternately exposed to aeration and nutrients from the culture medium. The gentle vertical movement of the medium creates a dynamic interface between air and liquid at the cell surface, ensuring an environment with minimal shear stress, high aeration and nutrient availability, no foaming, and unrestricted oxygen supply. This efficient transfer of nutrients and oxygen effectively replicates the alternating exposure to nutrients and oxygen that cells experience in their native environment.

## Macroporous Carriers

Anchorage-dependent cells require a surface for growth, traditionally provided by T-flasks or roller bottles. However, scaling up these methods is challenging, leading to the development of macroporous carriers. These matrices support the attachment, growth, and proliferation of adherent cell lines, including animal, mammalian, and insect cells. Commonly used in packed-bed and Tide Motion® bioreactors, they offer scalable solutions for large-scale cell culture. With diverse physical and chemical properties, the choice of carrier depends on the specific anchorage requirements of the cultured cells.



### BioNOC II®

BioNOC II® is a macroporous carrier that supports the growth of anchorage-dependent cells including **animal, mammalian, and insect cells** in either serum-containing or serum-free culture media.

### Features

- Made of 100% PET
- Complies to <85>, <87>, <881>, ISO 10993-5: 2009
- High porosity (50 - 200 µm between fiber), and biocompatible to coating reagents for better attachment
- High surface area\* up to 15,000 cm<sup>2</sup> per 0.1 L packed bed volume for Vero cells
- Enhanced hydrophilicity, coating factors can be added



### BioMESH®

BioMESH® is a low-lint macrocarrier designed specifically for cell harvesting applications. It is suitable for use with **skin fibroblasts, cardiac fibroblasts, multipotent/pluripotent stem cells, epithelial cells, and chondrocytes**.

### Features

- Made of a combination of polypropylene (PP) netting and PET mesh
- Complies to USP <83> <87> <788>, USP Class VI, ISO 10993-5: 2009
- High porosity (200 µm space between PET fiber; 670 µm space between PP fiber)
- High surface area\* for cell growth: up to 10,000 cm<sup>2</sup> per 0.1 L packed bed volume for MSCs
- Enhanced hydrophilicity, coating factors can be added


\*Note: Surface area varies on the cell line used



### Your Own Scaffold or Bioedible Scaffold

Choose your preferred carrier system or opt for a bioedible scaffold designed for seamless integration.





## Key Benefits

### Highest Yield

MiniTide® employs advanced macroporous carriers, BioNOC II® and BioMESH®, to optimize cell culture applications. BioNOC II® offers a 100% PET, fibrous culture matrix, while BioMESH® utilizes a combination of polypropylene (PP) netting and PET mesh. These macrocarriers offer several advantages:

### High surface area:

The macroporous design provides cells with a significantly increased surface area for optimal attachment and growth.

### Biomimetic environment:

The 3D structure closely mimics a cell's natural in vivo environment, promoting healthy cell behavior.

### Efficient mass transfer:

The matrix vessel facilitates efficient exchange of oxygen and nutrients, ensuring optimal cell proliferation and high bioproduct yields.

## Affordable Cost

MiniTide® scaled down bioreactor also offers significant advantages for research applications utilizing precious starting materials, high-cost culture media, or expensive analytical techniques. This translates to a substantial reduction in the overall cost per experiment.

## Linearly Scalable Quality

Esco Adherent Cell Tide Motion® platform offers a groundbreaking solution for adherent cell cultures. This innovative, world-first packed bed bioreactor system boasts full linear scalability. Starting with research-scale units

designed for standard CO<sub>2</sub> incubators, the platform seamlessly expands to production volumes while maintaining the core Tide Motion® principle. This ensures consistent and reliable cell growth throughout the development and manufacturing process.

Cells	Estimated Cell Number per MiniTide Vessel	Carriers
MSC	$1.2 \times 10^6$	BioNOC II®**
	$1.3 - 2.8 \times 10^7$	BioMESH®***
IPSC	$1.2 - 3.1 \times 10^6$	BioNOC II®
Fibroblast	$1.2 \times 10^7$	BioNOC II®
Vero	$6.2 \times 10^6$	BioNOC II®
HEK293	$8.8 \times 10^7$	BioNOC II®
MDCK	$3.5 \times 10^7$	BioNOC II®
BHK-21	$1.5 \times 10^8$	BioNOC II®
Sf-9	$1.7 \times 10^8$	BioNOC II®
Sf-21	$1.2 \times 10^7$	BioNOC II®
HuH-7	$7.1 \times 10^7$	BioNOC II®

\*The density of MSCs varies depending on the tissue of origin, the donor's age and characteristics, and the composition of the culture medium.

\*\*BioNoc II® is optimally suited for the culture of MSCs when the goal is to obtain secretomes or small exosomes as the final products.

\*\*\*BioMESH® is optimally suited for the cultivation of MSCs when the primary goal is cell harvesting as the final product. However, in cases where it is essential to maintain the cells' morphology, particularly if they have been extensively characterized in a 2D environment, BioMESH® may be more appropriate for harvesting small exosomes.

*The choice of carrier type will ultimately depend on the specific requirements of the client's final experiment.*

# Environmental Control

For both normoxic and hypoxic cell cultures, temperature is carefully regulated using Esco CO<sub>2</sub> incubators or compatible third-party models. To maintain aseptic conditions, these incubators should ideally have HEPA or ULPA filters. Additionally, a cooling function is recommended to counteract the heat generated by MiniTide®




## Top Holding Time

**Top holding time** refers to the phase in the Tide Motion® cycle where macroporous carriers are fully submerged in the culture medium, allowing for efficient nutrient exchange. During this stage, essential nutrients and growth factors diffuse into the carriers, while waste products produced by the cells are removed. This process ensures a continuous supply of nutrients to support cell growth and metabolism, optimizing the culture conditions for anchorage-dependent cells.


## Bottom Holding Time

**Bottom holding time** refers to the phase in the Tide Motion® cycle where the macroporous carriers are exposed and transition into the oxygenation stage. During this phase, the oxygen or suppressed oxygen diffuse from the headspace within the vessel to the surface of the macroporous carriers.

*Approximately 20% of the culture medium remains in contact with the macroporous carriers throughout the phase. This facilitates the passive diffusion of oxygen, or limited oxygen in the case of hypoxic cultures, from the surrounding air into the cell layer adjacent to the carriers, ultimately reaching the 3D cell culture within the macroporous matrix.*



Application/ Macrocarriers	EV/Exosomes	Cell Therapy <i>(when harvesting the cells as final product)</i>	Intracellular Virus <i>(when harvesting the cells containing the virus)</i>	Secreted Bioproducts <i>(Virus, Viral Vectors, Proteins)</i>
BioNOC II®	+++	++	++	+++
BioMESH®	+++	+++	+++	++




### BioNOC II® / BioMESH® (Packed Bed)

Cells remain entrapped in the carriers, simplifying media replacement and product harvesting.

### Cable Complete Set

Connects the base frame to the main control box.



### MiniTide® Culture Stage

Stainless steel 304 that can hold up to 4 MiniTide® vessels.

### HDMI Box

MiniTide® HMI utilizes Siemens HMI KTP700/PLC 1214C, which allows for independent control of the experiment conditions.



## Polystyrene


Polystyrene stands out as a preferred choice for cell culture flasks. This selection is driven by its key properties: superior optical clarity, ensuring unobstructed observation, exceptional durability for reliable performance, and minimal cytotoxicity.

## Polypropylene

Polypropylene exhibits exceptional resistance to moisture, a variety of acids, and alkalis. This material offers a robust combination of high impact strength and fatigue resistance, creating a favorable environment for cell growth and proliferation.

## Working Volume


The MiniTide® system employs a versatile cell culture vessel with a maximum capacity of 60 mL, but a working volume of 35 mL is recommended.



<b>Working Volume</b>	35 mL	
<b>Pre-Packed Carrier</b>	BioNOC II®	0.17 ± 0.005 g
		25 pcs
<b>Pre-Packed Carrier</b>		0.5 ± 0.05 g
	BioMESH®	Initially supplied as 10 BioMESH® sheets (24 × 9 mm). These sheets are bisected to yield a total of 20 BioMESH® sheets (12 × 9 mm).
<b>Material</b>	Cap	Polypropylene
	Bottle	Polystyrene

## MiniTide® HMI Main Screen

MiniTide-4 features independent control, enabling researchers to design diverse experiments. This allows for the exploration of various parameters, including independent manipulation of top holding (nutrient exchange) and bottom holding (aeration phase) durations, as well as the investigation of the impact of different media types.



HMI Main Screen	Displays the main screen
Trend	Displays the trend graph
Data	Displays the data list
Alarm	Displays the alarm list
Event	Displays the operation history
Setup	Displays the programming settings page and the manual operation page
Recipe	Displays recipe operation page
System	Displays the system configuration settings page
Log In	Displays user login screen
Log Out	Logging out the current user
Safety Remove Disk	Ensuring the files are properly saved
Safety Power Off	Function key to confirm closing the background program



# MiniTide® Product Specification

Esco Tide Motion® system bioreactors are a leading technology platform for global commercial production of human and animal vaccines. Additionally, Esco Tide Motion® system bioreactors play an essential role in clinical applications, supporting the development of adherent cell therapies and their derivatives. These derivatives include exosomes, extracellular vesicles, cell-derived nanovesicles, collagen for medical and cultivated meat applications, proteins for diagnostics and therapeutics, growth factors, and biomaterials like cultivated leather, all with applications in the biomedical field. MiniTide® is for research use or further manufacturing only and cGMP/GMP usage is subject to the site Quality Manufacturing Systems to be in place.

*Note: MiniTide® needs to be placed in a cooled CO<sub>2</sub> incubator to offset the heat it produces.*

<b>Weight</b>	HDMI Box**	6.5 kg
	Main control box + 1 MiniTide® Module unit	11 kg
	Main control box + 2 MiniTide® Module unit	13.7 kg
	Main control box + 3 MiniTide® Module unit	16.7 kg
	Main control box + 4 MiniTide® Module unit	20 kg
<b>Required space (W × H × D)</b>	HDMI Box**	290 × 360 × 290 mm
	Main control box + 1 MiniTide® Module unit	210 × 270 × 460 mm
	Main control box + 2 MiniTide® Module unit	295 × 270 × 460 mm
	Main control box + 3 MiniTide® Module unit	380 × 270 × 460 mm
	Main control box + 4 MiniTide® Module unit	465 × 270 × 460 mm
<b>Power supply</b>	100-230 VAC, 50/60 Hz	
<b>Housing</b>	304 stainless steel	
<b>Tide Motion® rate range</b>	15 ~ 75 mL/min	
<b>Accuracy</b>	1.5 mL/min <sup>1</sup>	
<b>Holding Time</b>	0 ~ 9999 s	
<b>Other properties</b>	Data trend chart Data sheet Alarm list Recipe management User authority management Export in PDF and CSV format	

\* The user needs to ensure that the third-party CO<sub>2</sub> incubator with a cooling function can withstand the weight

\*\* The HMI box is placed outside of the CO<sub>2</sub> incubator

Item Code	Model Code	Description
<b>Main Unit</b>		
2231110	MIT-1M10	1 VESSEL, MODULE, 110-230V, 50/60HZ, - M for Modular
2231111	MIT-2M10	2 VESSELS, MODULE, 110-230V, 50/60HZ
2231112	MIT-3M10	3 VESSELS, MODULE, 110-230V, 50/60HZ
2231113	MIT-4M10	4 VESSELS, MODULE, 110-230V, 50/60HZ
2231097	MINITIDE, MIT-16A-8	16 VESSELS, AUTO, 230V, 50/60HZ

*Note: Existing MiniTide® from stages 1 through 4 are capable of being retrofitted as per user's requirement.*

## MiniTide® Product Specification

Consumables		
MiniTide Vessel (MTMV)		
1400309	MINITIDE VESSEL-E, MTMV-E	4 EMPTY VESSELS *Matrix depends on client's requirement/s
1400310	MINITIDE-BIONOC II, MTMV-BII, 10 PCS	For virus and applications which are sensitive to pH 10 pieces of BioNOC II Matrix:Media Ratio is 1:20 4 Vessel/Pack
1400311	MINITIDE-BIONOC II, MTMV-BII, 25 PCS	25 pieces of BioNOC II Matrix:Media Ratio is 1:7 4 Vessel/Pack
1400312	MINITIDE-BIOMESH 50 MICRON, MTMV-BM50	MiniTide-BioMESH 50 micron Matrix:Media Ratio is 1:7 4 Vessel/Pack
1400313	MINITIDE-BIOMESH 200 MICRON, MTMV-BM200	MiniTide-BioMESH 200 micron Matrix:Media Ratio is 1:7 4 Vessel/Pack
MiniTide Base Frame		
5012040	BASE FRAME, MINITIDE-2	SS CARCASS ONLY
5012041	BASE FRAME, MINITIDE-3	SS CARCASS ONLY
5012042	BASE FRAME, MINITIDE-4	SS CARCASS ONLY

## North America

### USA

#### Esco Technologies, Inc.

##### Life Sciences/Medical

903 Sheehy Drive, Suite F, Horsham, PA 19044, USA  
Toll-Free USA and Canada:  
1-877-479-3726  
Tel: +1 215 441 9661  
Fax: 484 698 7757  
Email: eti.admin@escolifesciences.com  
Website: www.escolifesciences.com

##### Healthcare Division/Factory

2512 Metroprolitan Dr, Suite 120-B  
Feasterville-Treove, PA 19053-6738  
Tel: +1 215 322 2155  
Email: eti.pharma@escolifesciences.com  
Website: www.esco-pharma.com

## Europe

### DENMARK

#### Esco Medical ApS

Kringelleg 10, DK-8250 Egaa, Denmark  
Tel: +45 5397 3067  
Email: medical@escolifesciences.com

### UK

#### Esco GB Ltd

Unit 2 R-Evolution @ Gateway 36,  
Kestrel Way, Barnsley S70 5SZ  
Tel: +44 (0) 1226 360 799 (Pharma)  
+44 (0) 1226 361 529 (Lab)  
Email: egb.info@escolifesciences.com  
Websites: www.esco-pharma.com  
www.escolifesciences.com

### GERMANY

#### Esco Lifesciences GmbH

Straßheimer Straße 17  
61169 Friedberg, Germany  
Tel: +49 6031 6873447  
Email: ede.info@escolifesciences.com  
Website: www.esco-global.de

### RUSSIA

#### Esco Russia

Building 4, 8, Novovladikinskii proezd  
Moscow, Russian Federation 127106  
Email: esco.russia@escolifesciences.com  
Website: escolifesciences.ru

### LITHUANIA

#### Esco Medical Technologies UAB

Gamybos g. 2, Ramuciai, Kauno r., 54468  
Lithuania  
Email: medical@escolifesciences.com  
support-medical@escolifesciences.com  
Website: esco-medical.com

## Asia Pacific

### BANGLADESH

#### Esco Lifesciences (Bangladesh) Pvt. Ltd.

H # 662 (6th Floor), Apt. # A-6, R # 9  
Mirpur DOHS, Mirpur, Dhaka-1216,  
Bangladesh  
Tel: +8801958541100  
Email: rafia.oshin@escolifesciences.com  
ebd.customerservice1@escolifesciences.com

### CHINA

#### Esco Shanghai Trading Co. Ltd.

Email: mail@escolifesciences.cn  
Website: www.escolifesciences.cn

##### Beijing

Rm. 502, Tower A, Times Fortune,  
Sanyuanqiao, Chaoyang District, Beijing,  
China 100028  
Phone: +86 (10) 5867 7868 / 5867 7823  
Fax: +86 (10) 5867 9244  
Email: dora.miao@escolifesciences.com

##### Guangzhou

Room 802, Yuanyang Building,  
20 Huacheng Avenue, Tianhe District,  
Guangzhou  
Tel: +86 (20) 3837 3621

##### Qingdao

Room 1809, Cyberport Flagship Building,  
No. 40 Hong Kong Middle Road, Shinan  
District, Qingdao,  
Shandong Province, China  
Email: ailey.wang@escolifesciences.com

##### Shanghai

4-5th Floor, Building 16, No.388  
Shengrong Road, Pudong New District,  
Shanghai, China  
Tel: +86 (21) 5065 3007

##### Hong Kong (China)

#### Esco Technologies (Hong Kong) Limited

Unit 06-07, 27/F, New Treasure Centre,  
10 Ng Fong Street, San Po Kong,  
Kowloon, Hong Kong  
Tel: +852 3628 3986  
Fax: +852 3186 2821  
Email: mail@escolifesciences.hk  
Website: http://escolifesciences.hk

##### Taiwan (China)

**Esco Lifesciences Taiwan Co., Ltd.**  
No. 122-2, Wuqing Rd., Dayuan Dist.,  
Taoyuan City 337014, Taiwan  
Tel: +886 33818837  
Fax: +886 33817135  
Email: mail@escolifesciences.tw  
Website: http://escolifesciences.tw

### INDIA

#### ELG Lifesciences India Private Limited

2nd Floor, No.7/3, KVD Towers, Old  
Madras Road, Kolar Road, Bengaluru,  
Karnataka, 560038  
Tel: +91 8971532303  
Email: india@escolifesciences.com

### INDONESIA

#### PT Esco Utama

Taman Tekno BSD Blok H8 No.1 BSD  
Sektor XI, Serpong, Tangerang Selatan  
15314  
Tel: +6221 30 446 565  
Fax: +6221 30 446 566  
Email: eucs2@escolifesciences.com  
Website: http://escolifesciences.co.id

### MALAYSIA

#### Esco Micro (M) Sdn. Bhd.

**Technical Support Hotline:**  
1 300 88 ESCO (3726)

##### Selangor

No. 15, Jalan Sungai Buloh 27/101A,  
Persiaran Klang, HICOM 27 Industrial  
Park,  
40400 Shah Alam, Selangor, Malaysia  
Tel: +603 5103 8480  
Fax: +603 5614 3385  
Email: malaysia@escolifesciences.com  
Website: https://www.escolifesciences.com

### Melaka

No 109 & 109-1, Jalan Zapin 1,  
Taman Paya Rumput Utama,  
75460, Melaka, Malaysia  
Tel: 010-2638007  
Website: https://www.escolifesciences.com

### MYANMAR

#### Esco Lifesciences Co.,Ltd

No. 376 (E), Shwe Hnin Si Road, Say Ou  
Yin Street, 5 Ward, Mayangone Township,  
Yangon, Myanmar  
Tel: +959 970 615 645  
+959 961 456 426  
Email: ingyin.haymarn@escolifesciences.com  
eaint.thaemon@escolifesciences.com

### PHILIPPINES

#### Esco Philippines, Inc.

##### Manila

Unit 707E, 7th floor, East Tower Four  
E-com  
Bldg, Block 22 Seaside Cor Diokno Ave.  
MOA Complex, Pasay City 1300  
Direct Line: +63 (02) 8828 1527  
Trunkline: +63 (02) 8478 0384  
Email: philippines@escolifesciences.com  
Website: https://www.escolifesciences.com

##### Cebu

Unit 2K, Baseline Residences, Juana  
Osmena St. Brgy. Capitol Site, Cebu City  
6000  
Tel: +63 (32) 232 5826

##### Davao

71-C Lupo Diaz St. cor Nicasio Torres St.,  
Bo. Obrero, Davao City 8000  
Tel: +63 (082) 234-0997

### SINGAPORE

#### Esco Micro Pte Ltd

19 Changi South Street 1,  
Singapore 486779  
Tel: +65 65420833  
Fax: +65 65426920  
Email: mail@escolifesciences.com  
Website: https://www.escolifesciences.com

### SOUTH KOREA

#### Esco Korea Micro Pte Ltd

2F, 845, Geumo-ro, Gwangmyeong-si,  
Gyeonggi-do, Republic of Korea (14297)  
Tel: +82 2 830 0482  
Fax: +82 2 830 0491  
Email: info@escoglobal.co.kr  
Website: http://www.esco-global.co.kr

### THAILAND

#### Esco Lifesciences (Thailand) Co.,Ltd.

8/3 Soi Rimthangduan 2, Sukhumvit Road,  
Bangchak,  
Phrakhanong, Bangkok 10260  
Tel: 02 082 2029  
Fax: 02 117 3746  
Email: cst.th.sales1@escolifesciences.com  
Official Line: @escothailand  
Website: http://escolifesciences.co.th/

### VIETNAM

#### Esco Vietnam Company Ltd.,

Email: vietnam@escolifesciences.com  
Website: vn.esco-global.com

##### Hanoi

Floor 2, Dream Home Center Building,  
282 Nguyen Huy Tuong street, Thanh Xuan  
Trung Ward, Thanh Xuan District, Hanoi City,  
Vietnam  
Tel: +84 24 6269 1460  
+84 24 6269 1461

### Ho Chi Minh

No. 14 - Lot G3, Street No. 53,  
Tan Quy Dong Settlements,  
Tan Phong Ward, District 7, HCMC, Vietnam  
Tel: +84 28 3776 0363

## Middle East

### DUBAI

#### Esco Lifesciences Trading LLC

Showroom 2 SW2 Gate 1 Jadaf 1 Dubai, UAE  
P.O Box: 182099  
Tel.: +971 (4) 5461137  
Email: gildence.casipit@escolifesciences.com  
abdalmageed.shalaby@escolifesciences.com

## Africa

### SOUTH AFRICA

#### Esco Technologies Pty Ltd

##### Centurion

Unit 2 Landmark Park,  
17 Landmarks Avenue,  
Kosmosdal Extension 11,  
Centurion 0157 South Africa  
Tel: +27 (0) 11 314 3184  
Mobile: +27 (0) 82 853 9655  
Email: sonja.strydom@escolifesciences.com  
Website: https://www.escolifesciences.com

##### Cape Town

Unit 2, 7 Edison Way, Century City,  
Cape Town, 7441  
Tel: +27 (0) 11 314 3184

## Oceania

### AUSTRALIA

#### Esco LifeSciences (Australia) Pty Ltd

Suite 7, Level 7, 330 Collins Street,  
Melbourne,  
Victoria 3000, Australia



**ESCO LIFESCIENCES GROUP NETWORK**  
*42 Locations in 21 Countries All Over the World*



- Global Offices
- Distributors
- Factories
- R&D Centers
- Regional Distribution Centers



- Air Shower
- Aseptic Containment Isolator (ACTI)
- Ceiling Laminar Airflow Units
- Cleanroom Transfer Hatch
- Containment Barrier Isolator (CBI)
- Downflow Booth (DFB)
- Dynamic Floor Laminar Hatch
- Dynamic Pass Box
- Evidence Drying Cabinet
- Garment Storage Cabinet
- General Processing Platform Isolator (GPPI)
- Laminar Flow Horizontal Trolley
- Laminar Flow Straddle Units, Single and Double
- Laminar Flow Vertical Trolley
- Pass Box
- Soft Wall Cleanroom
- Sputum Booth
- Ventilated Balance Enclosure (VBE)
- Weighing and Dispensing Containment Isolator (WDCI)

Since 1978, Esco has emerged as a leader in the development of controlled environment, laboratory and pharmaceutical equipment solutions. Products sold in more than 100 countries include biological safety cabinets, fume hoods, ductless fume hoods, laminar flow clean benches, animal containment workstations, cytotoxic cabinets, hospital pharmacy isolators, and PCR cabinets and instrumentation. With the most extensive product line in the industry, Esco has passed more tests, in more languages, for more certifications, throughout more countries than any biosafety cabinet manufacturer in the world. Esco remains dedicated to delivering innovative solutions for the clinical, life science, research and industrial laboratory community.



**Esco Healthcare**  
 19 Changi South Street 1 Singapore 486779  
 Tel: +65 65420833  
 Email: mail@vaccixcell.com



**Esco Lifesciences Offices:** Bangladesh | China | Denmark | Germany | Hong Kong | Indonesia | Italy | Lithuania | Malaysia | Myanmar | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam

Mintide Brochure\_A4\_v8\_Blue\_03262025  
 Esco Lifesciences Group reserves the right to alter its products and specifications without notice. All trademarks and logos in this material are the property of Esco Lifesciences Group and the respective companies.