



Standalone
CelCradle X[®]



CCX[®]-SAH

CelCradle X[®]

Closed Automated Single-use cGMP/
GMP 13485 Bioreactor with Closed
Single-use Harvester

CelCradle X[®]

Closed Automated Single-use cGMP/GMP 13485 Bioreactor with Closed Single-use Harvester

Features

- Complete standalone bioreactor running on Siemens HMI/PLC with single-use bottles, pH, and DO sensors
- Can be used as a standalone bioactor or integrated with an Esco Laminar Airflow Cabinet, Biosafety Cabinet, or Cell Processing Isolator
- Equipped with H₂O₂ biodecontamination system
- Produced in an ISO 9001/14001/13485 factory
- Automated seeding and control of pH/DO and glucose (optional)
- Single use automated harvesting system for autologous cell therapy, process development, and seed train for TideXcell[®]
- Able to use autoclavable pH/DO sensors for process development and switch *in situ* for single use pH/DO



CCX Isolator Integration:
designed as dock/undock bioreactor module to cell processing isolator



Standalone CCX on Support Stand



Commercial cGMP/GMP Solutions

- Vaccine Production
- Monoclonal Antibodies or Recombinant Proteins for therapeutics/*in vitro*/*in vivo* diagnostics
- Drug Discovery
- Conditioned media production for cosmeceuticals
- Adherent cell master cell bank (MCB) or working cell bank (WCB) generation
- Intracellular, Secreted Viral Vector or Virus Process Development
- Proteome Research
- PK/PD Testing
- Virus MVB production

Enabling research/clinical solutions

- Adherent cell seed train for TideXcell®
- Stem cell or exosome-autologous therapy (bone regeneration using bioscaffolds)
- Autologous adherent cell therapy (e.g. hair loss, face recontouring, etc.)
- Technology transfer from multi-stack or hyperstack bioreactors for small scale production of vaccines/virus/viral vectors
- Monoclonal antibodies or proteins
- Biochemistry/Stemistry



Options:



- Minus 8°C from ambient to 40°C
- Aquaculture vaccines (as low as 16 °C)
- Low temperature 25-28°C with 2.5% CO₂



- Mini fridge for final secreted product collection



- Hollow fiber for protein concentration and clarification



- Aseptic sampling



- Improve RH up to 95%



- Orbital shaker for insect or other cell seed train

NOMENCLATURE					
Model	Bioprocessing Mode	Bioprocessing Mode/ Bottles	Mixing Bottles	Packed Bed Volume	Carrier Type
CelCradle X®	SA-SUGH (Semi-Automated Single Use Glass Hybrid) Single use bottle with single pump, glass mixer, autoclavable pH/DO sensors	500AP Single use. Batch with Re-circulation 2L (0.5L cannot be used)	Single use 2L bottle, single use pH & DO sensor	Empty	Client's own or pack lower volume BioNOC™ II/BioMesh®
	SA-SUB (Semi-Automated Single Use CC Bottle with peristaltic pump) After cells have reached confluency, the CC bottle will be disconnected and placed into a separate CelCradle Harvester equipped with in-line, purification, filtration, and filling	500-AP Single use. Perfusion with Feed/Harvest of Biologic 0.5L. Mixer (2L cannot be used), 2L Feed / Harvest Bag	Single use 0.5L, single use pH & DO sensor	0.1 liter	BioNOC™ II/BioMesh®
	FA-SUB (Fully Automated Single Use from cryovial to filled product) <i>*under development</i>		Glass 2L autoclavable pH & DO sensor	0.2 liter <i>*not available for beta testing</i>	BioNOC™ II/BioMesh®
			Glass 0.5L autoclavable pH & DO sensor	0.3 liter <i>*not available for beta testing</i>	BioNOC™ II/BioMesh®

**Specifications subject to change.*

CelCradle X[®]-SAH

The CelCradle X[®] Semi-Automated Harvesting System (CCX-SAH) is used to harvest whole cells from the CelCradle bottle through automation by setting up various parameters.

This semi-automated cell harvesting of adherent CelCradle cultures includes but not limited to: Insect Cells, Mammalian Cells, Stem Cells.

The mounted fixed-speed peristaltic pumps provides for convenient liquid addition/removal and cell suspension harvest during the culture's harvesting process.

Applications

- Autologous cell therapy
- Mammalian cell harvest
- Biomass harvesting for seed cell / seed train
- Biomass harvesting for 3D bioprinting
- Biomass harvesting for cellular agriculture/lab grown meat
- Harvesting automation for secreted or non-secreted products.
- Harvesting for intracellular viruses

Features

- Mechanical dissociation of adherent cells from Tide Motion macroporous carriers
- Automated harvesting run can be setup from washing to harvesting (secreted or non-secreted products) in a closed, automated manner.
- Suitable for cell harvest from CelCradle X[®] System
- Support the cell seed source for TideXcell[®]-002 up to TideXcell[®]-100 systems.
- Harvest viable animal cells within 1 to 3 hours (harvest time depends on harvest cycles).
- Able to achieve more than 90% harvesting efficiency

Note: Cell harvest efficiency depends on the cell line and dissociation enzyme used.

Parts of CCX®-SAH



CelCradle X® Semi-Automated Harvesting System (CCX-SA)

Exterior dimension (WxDxH)	SAH: 928 x 380 x 822 mm (36.53 x 14.96 x 32.36 in) Incubator: 543 x 640 x 759 mm (21.38 x 25.20 x 29.88 in)
Weight	SAH: 50 kg (110.23 lb) Incubator: 85 kg (187.39 lb) Whole Unit: 135 kg (297.62 lb)
Power	Single phase, AC, 220 V, 50 Hz
Motor	Shaker motor
Max Speed range	0-300 rpm
Tubing	ID 9.6 mm(3/8")*WT 3.2 mm(1/8")
Capacity	1 CelCradle bottle
Consumable CelCradle Bottle Volume	500 ml
CelCradle Material	PP
Incubator temperature range	20-38°C
Touchscreen Size	Control panel: 7" color Incubator HMI: 7" color

*Specifications subject to change.

ESCO LIFESCIENCES GROUP GLOBAL OFFICES

North America

USA

Esco Technologies, Inc.

Life Sciences/Medical

903 Sheehy Drive, Suite F, Horsham,
PA 19044, USA

Toll-Free USA and Canada:

1-877-479-3726

Tel: 215-441-9661

Fax: 484-698-7757

Email: eti.admin@escoglobal.com

Website: <http://escolifesciences.us>

Healthcare Division/Factory

2512 Metropolitan Dr, Suite 120-B
Feasterville-Treose, PA 19053-6738

Tel: +1 215-322-2155

Email: eti.pharma@escoglobal.com

Website: www.escopharma.com

Europe

DENMARK

Esco Medical ApS

Kringelled 10, DK-8250 Egaa, Denmark

Tel: +45 5397 3067

Email: medical@escoglobal.com

UK

Esco GB Ltd

Unit 2 R-Evolution @ Gateway 36,

Kestrel Way, Barnsley S70 5SZ

Tel: +44 (0) 1226 360 799 (Pharma)

+44 (0) 1226 361529 (Lab)

Email: egb.info@escoglobal.com

Websites: www.escopharma.com

www.escolifesciences.co.uk

GERMANY

Esco Lifesciences GmbH

Straßheimer Straße 17

61169 Friedberg, Germany

Tel: +49 6031 6873447

Email: mail@escoglobal.com

Website www.escoglobal.de

ITALY

Esco Technologies S.r.l.

Operative and Legal office:

Via Poggio Moiano, 23,

00199, Rome, Italy

Tel: +39 06 69271849

Email: ilenia.nardolillo@escoglobal.com

Website: escolifesciences.it

RUSSIA

Esco Russia

Building 4, 8, Novovladikinskii proezd

Moscow, Russian Federation 127106

Email: esco.russia@escoglobal.com

Website: escolifesciences.ru

LITHUANIA

Esco Medical Technologies UAB

Draugystes 19

51230 Kaunas, Lithuania

Website: escomedicalgroup.com

Email: medical@escoglobal.com

support-medical@escoglobal.com

Asia Pacific

BANGLADESH

Esco Lifesciences (Bangladesh) Pvt. Ltd.

H # 662 (6th Floor), Apt. # A-6, R # 9

Mirpur DOHS, Mirpur, Dhaka-1216,

Bangladesh

Tel: +8801907-700777

Mobile: +88 01911 629400

Email: ebd.customerservice1@escoglobal.com

CHINA

Esco Shanghai Trading Co. Ltd.

Email: mail@escolifesciences.cn

Website: www.escolifesciences.cn

Beijing

Rm.502, Tower A, Times Fortune,

Sanyuanqiao

Chaoyang District, Beijing

Tel: +86 (10) 5867 7868

Fax: +86 (10) 5867 9244

Shanghai

Room1211, Jiahe International Building,

No.1, Lane 66, Huayuan Road,

Hongkou District, Shanghai

Tel: +86 (21) 6095 1955

Guangzhou

Room 901, West Block, No. 21,

Huali Road, Tianhe District, Guangzhou,

China 510623

Tel: +86 (20) 3837 3621

Chengdu

Rm. 1-406, Building No. 2, Shangding

International Building,

No.27 4th Block, South Renmin Road,

Chengdu City,

Sichuan Province

Tel: +86 (28) 8553 6219

Qingdao

Room 23C, Fucai Building, No.135 of

Yan'an Road,

Shinan District, Qingdao City,

Shandong Province

Tel: +86 (532) 8163 9321

HONGKONG

Esco Technologies (Hong Kong) Limited

Unit 904, Laurels Industrial Centre,

No.32 Tai Yau Street, San Po Kong,

Kowloon, Hong Kong

Tel : +852 3628 3986

Fax : +852 3186 2821

Email: mail@escolifesciences.hk

Website: <http://escolifesciences.hk>

INDONESIA

PT Esco Utama

Taman Tekno BSD Blok H6 no 10,

Serpong, Tangerang 15314

Tel: +6221 29 666 565 / 301

Fax: +6221 29 666 302

Email: eucs2@escoglobal.com

Website: <http://escolifesciences.co.id>

MALAYSIA

Esco Micro (M) Sdn. Bhd.

Technical Support Hotline:
1 300 88 ESCO (3726)

Selangor

No. 15, Jalan Sungai Buloh 27/101A,
Persiaran Klang, HICOM 27 Industrial
Park,
40400 Shah Alam, Selangor, Malaysia
Tel: +603 5103 8480
Fax: +603 5614 3385
Email: malaysia@escoglobal.com
Website: http://my.escoglobal.com

Melaka

No. 13-1, Jalan Melaka Raya 32,
Taman Melaka Raya 1,
75000 Melaka, Malaysia
Tel: +606 284 0007
Fax: +606 283 0288
Website: http://my.escoglobal.com

MYANMAR

Esco Lifesciences Co.,Ltd

717(B), Bayintnaung Road,
Insein Township, Yangon, Myanmar
Tel: +959 423 488 984
+959 963 520 535
Email: ingyin.haymarn@escoglobal.com
csts4@escoglobal.com

PHILIPPINES

Esco Philippines, Inc.

Manila

Unit 710E, 7th floor,
East Tower- Four Ecom
Bldg, Block 22 Seaside Cor Diokno Ave.
MOA Complex, Pasay City 1300
Direct Line: +63 (02) 828 1527
Trunkline: +63 (02) 478 0384
Email: philippines@escoglobal.com
Website: http://escoglobal.com.ph

Cebu

Blk 3 Lot 13 Mutual Homes Subd.
Pajac, Lapu-Lapu City, Cebu 6015
Tel : +63 (32) 232 5826

Davao

Unit 104 & 105, G/F Metro Lifestyle
Building
F. Torres St. Cor. E Jacinto Ext., Davao City
8000
Tel : +63 (82) 333 4132

SINGAPORE

Esco Micro Pte Ltd

21 Changi South Street 1
Singapore 486777
Tel: +65 65420833
Fax: +65 65426920
Email: mail@escoglobal.com
Website: www.escoglobal.com

South Korea

Esco Korea Micro Pte Ltd

#206, Byucksan Digital Valley 8-Cha, 70,
Gyeongin-ro 71-gil, Yeongdeungpo-gu,
Seoul,
Republic of Korea
Tel: +82-2-830-0482
Fax: +82-2-830-0491
Email: info@escoglobal.co.kr
Website: http://www.escoglobal.co.kr

TAIWAN

Esco Lifesciences Taiwan

No. 112, Sec. 2, Qingfeng Rd., Zongli Dist.,
Taoyuan City 320, Taiwan
Tel: +886 34533220
Fax: +886 34532770
Email: mail@escolifesciences.tw
Website: http://escolifesciences.tw

THAILAND

Esco Lifesciences (Thailand) Co.,Ltd.

8/3 Soi Rimthangduan 2, Sukhumvit Road,
Bangchak, Phrakanong , Bangkok 10260
Tel: 062 813 3331
Fax: 02 117 3746
Email: csts.th.sales1@escoglobal.com
Website: http://escolifesciences.co.th/

VIETNAM

Esco Vietnam Company Ltd.,

Email: vietnam@escoglobal.com
Website: vn.escoglobal.com

Hanoi

7th Floor, No. 8, Alley 15,
Trung Kinh Street,
Trung Hoa Ward, Cau Giay District,
Hanoi, Vietnam
Tel: +84 24 6269 1460
+84 24 6269 1461

Ho Chi Minh

No. 14 - Lot G3, Street No. 53,
Tan Quy Dong Settlements,
Tan Phong Ward, District 7, HCMC,
Vietnam
Tel: +84 28 3776 0363

Middle East

DUBAI

Esco Lifesciences Trading LLC

Dubai
Esco Lifesciences Trading LLC
1529 Tamani Arts Offices, Al Asayel St.,
Business Bay, Dubai UAE
P.O Box: 182099
Tel: (04) 770 6674
Email: alvin.heah@escoglobal.com

Africa

SOUTH AFRICA

Esco Technologies Pty Ltd

Unit 2 Landmark Park,
17 Landmarks Avenue,
Kosmosdal Extension 11,
Centurion 0157 South Africa
Tel : +27 (0) 11 314 3184
Mobile: +27 (0) 82 853 9655
Email : sonja.strydom@escoglobal.com
Website: http://za.escoglobal.com

ESCO[®]
Healthcare



ESCO
VACCIXCELL
ADHERENT BIOPROCESSING SPECIALIST

21 Changi South Street 1 • Singapore 486777
Tel +65 6542 0833 • Fax +65 6542 6920 • mail@vaccixcell.com
www.escovaccixcell.com

CDMO Services:
Esco Aster Pte Ltd
#02-09, 67 Ayer Rajah Crescent, Singapore 139950
Tel +65 6251 9361 • mail@escoaster.com
www.escoaster.com